

Aegispak INNO Column INNO-P Column INNO-Core Column

High Performance Liquid Chromatography Column

Ghost Zero Column

Solvent Filter



YJ BIOCHROM

<http://www.yjbiochrom.com>

Product Group


Aegispak Columns

High purified silica gel surface was coated with silicone polymer by chemical vapor deposition (CVD), and the functional group was directly bonded from coated silica gel.

(Gas phase reaction synthesis method)

After the resin synthesis, End-capping was performed at high temperature and high pressure to remove the remaining silanol groups. This minimizes the effect of residual silanol groups to the extreme.


This resin is functional group directly bonded to a silicone polymer(-Si-R). There is a slight difference in the separation pattern from ordinary resin (-Si-O-Si-R).



INNO Columns

This resin was synthesized by directly bonding the functional group to the silanol group of refined high purity silica gel. (-Si-O-Si-R) (liquid phase reaction synthesis method)

End -capping was performed 2 times at a high temperature to remove the remaining silanol groups after the resin synthesis. Since all manufacturing processes are carried out in solution, it is possible to bulk production.



INNO-P Columns

This resin was synthesized by directly bonding the functional group to the silanol group of refined high purity silica gel with pore size of 300Å. (-Si-O-Si-R) (liquid phase reaction synthesis method)

End -capping was performed at a high temperature to remove the remaining silanol groups after the resin synthesis.



INNO-Core Columns

This column contains a non-porous solid core in porous silica gel. This column is capable of higher-speed analysis than common silica gel-based columns.



Ghost Zero Column

This Column is a filter that removes unknown peaks or ghost peaks generated in the mobile phase.

This Column is used by connecting between the pump and the auto-sampler.




| Classify | Aegispak Column | INNO Column | INNO-P Column | INNO-Core Column |
|---------------|-----------------|-------------|---------------|------------------|
| C18 | • | • | • | • |
| C18-F | • | | | |
| C18-FA | • | | | |
| C18-L | • | | | |
| C18-LA | • | | | |
| C18-SB | | • | | |
| C18-SSB | | • | | |
| C18-SBA | | • | | |
| C18-HC | | • | | |
| C18-NE | | • | | |
| C18-PE | | | | • |
| C8 | • | • | • | • |
| C8-SB | | • | | |
| C8-HC | | • | | |
| C8-NE | | • | | |
| C4 | | • | • | |
| C1 | | • | | |
| Phenyl | | • | | |
| Phenyl-Hexyl | | | | • |
| PFP | | | | • |
| NH2 | | • | | |
| CN | | • | | |
| Silica | | • | | |
| Diol | | • | | |
| SCX | | • | | |
| SCX-L | | • | | |
| SAX | | • | | |
| C18/SCX (1/1) | | • | | |
| C18/SCX (5/1) | | • | | |
| C18/SAX (1/1) | | • | | |
| C18/SAX (5/1) | | • | | |

| Classify | HPLC | UPLC |
|-------------------|------|------|
| Ghost Zero Column | • | • |




Evaluation of C18 Resin

■ Amitriptyline




- Amitriptyline is a strong basic compound
- Our columns have less tailing and good shape.

Aegispak C18-L




INNO Column C18



(HPLC Condition)

- **Column Size** : 4.6 X 150mm, 5 μ m
- **Flow Rate** : 1.0 ml/min
- **Temperature** : 35 $^{\circ}$ C
- **Mobile Phase** : Methanol / 20mM Potassium Phosphate (Mono + Dibasic), pH 7.0 = 8:2
- **Samples** : A) Uracil
B) Toluene
C) Ethylbenzene
D) Quinizarin
E) Amitriptyline

■ Phenylacetylacetone



- Phenylacetylacetone is a coordination compound.
- Our columns have less tailing and good shape.

Aegispak C18-L







INNO Column C18




(HPLC Condition)

- **Column Size** : 4.6 X 150mm, 5 μ m
- **Flow Rate** : 1.0ml/min
- **Detection** : PDA 254nm
- **Temperature** : 40 $^{\circ}$ C
- **Mobile Phase** : Methanol / DW = 5:5
- **Samples** : A) Uracil
B) Caffein
C) Phenol
D) 2-Ethyl Pyridine
E) Methyl Benzoate
F) Benzene
G) Toluene
H) Phenylacetylacetone
I) Naphthalene

Retention Time Character




INNO Column Silica 100Å



(HPLC Condition)

- **Column Size** : INNO Column Silica 4.6 x 250mm 5µm
- **Flow Rate** : 1.0 mL/min
- **Detection** : PDA 254 nm
- **Temperature** : 35°C
- **Mobile Phase** : n-Hexane / Ethyl alcohol = 98 / 2
- **Samples** : 1) Toluene : 8.8mg/ml
2) p-Dinitrobenzene : 2.2mg/ml


INNO Column NH2



(HPLC Condition)

- **Column Size** : INNO Column NH2 4.6 x 250mm 5µm
- **Flow Rate** : 1.0 mL/min
- **Detection** : RID
- **Temperature** : 35°C
- **Mobile Phase** : Acetonitrile / DW = 75 / 25
- **Samples** : 1) Fructose : 12 mg/mL
2) Glucose : 20 mg/mL
3) Saccharose : 20 mg/mL
4) Maltose : 30 mg/mL
in acetonitrile 50(Vol)%
- **Shipping Solvent** = 100% Acetonitrile

INNO Coulmn SCX / INNO Coulmn SCX-L




(HPLC Condition)

- **Column Size** : 4.6 X 250mm, 5µm
- **Flow Rate** : 1.0 mL/min
- **Detection** : PDA 254 nm
- **Temperature** : 35°C
- **Mobile Phase** : 0.2 M (NH4)H2PO4 (pH 3.5)
- **Samples** : 1) Uracil : 0.2mg/mL
2) Thymin : 0.3mg/mL
3) Cytosine : 0.6mg/mL




INNO Column SAX



(HPLC Condition)

- **Column Size** : INNO Column SAX
4.6 x 250mm, 5 μ m
- **Flow Rate** : 1.25mL/min
- **Detection** : UVD 254nm
- **Temperature** : 35°C
- **Mobile Phase** : Isooctane /Ethanol /DW
= 85 / 14.7 / 0.3
- **Samples** : 1) Toluene : 6.25mg/mL
2) Nitrobenzene : 0.63mg/mL
3) 2-Nitroaniline : 1.25mg/mL
4) 3-Nitroaniline : 1.25mg/mL
5) 4-Nitroaniline : 2.5mg/mL
- **Shipping Solvent** : 100% Ethanol



INNO Column C18/SCX



(HPLC Condition)

- **Column Size** : 4.6 x 250mm, 5 μ m
- **Flow Rate** : 1.0mL/min
- **Detection** : UVD 254nm
- **Temperature** : 35°C
- **Mobile Phase** : Acetonitrile /DW = 60 / 40
- **Samples** : 1) Uracil : 0.1mg/mL
2) Methyl benzoate : 2.2mg/mL
3) Toluene : 8.7mg/mL
4) Naphthalene : 0.9mg/mL
- **Shipping Solvent** : 60% Acetonitrile

INNO Column C18/SAX




(HPLC Condition)

- **Column Size** : 4.6 X 250mm, 5 μ m
- **Flow Rate** : 1.0 mL/min
- **Detection** : UVD 254nm
- **Temperature** : 35°C
- **Mobile Phase** : Acetonitrile /DW = 60 / 40
- **Samples** : 1) Uracil : 0.1 mg/mL 2) Methyl benzoate : 2.2 mg/mL 3) Toluene : 8.7 mg/mL
4) Naphthalene : 0.9 mg/mL
- **Shipping Solvent** : 60% Acetonitrile




INNO Column Diol



(HPLC Condition)

- **Column Size** : INNO Column Diol 4.6 X 250mm, 5 μ m
- **Flow Rate** : 1.0 mL/min
- **Detection** : PDA 254 nm
- **Temperature** : 35 $^{\circ}$ C
- **Mobile Phase** : n-Hexane / Ethyl alcohol = 85 / 15
- **Samples** : 1) 2-nitrotoluene 0.2mg/mL
2) Phenol : 1.0mg/mL
3) Resorcinol 2.5 mg/mL
4) hydroquinone 10mg/mL


INNO Column CN



(HPLC Condition)

- **Column Size** : INNO Column CN 4.6 X 250mm, 5 μ m
- **Flow Rate** : 1.0 mL/min
- **Detection** : PDA 254 nm
- **Temperature** : 35 $^{\circ}$ C
- **Mobile Phase** : Methanol / DW = 50 / 50
- **Samples** : 1) Uracil : 0.1mg/mL
2) Methyl benzoate : 2.2mg/mL
3) Toluene : 8.7mg/mL
4) Naphthalene : 0.9mg/mL

INNO Column Phenyl




(HPLC Condition)

- **Column Size** : INNO Column Phenyl 4.6 X 250mm, 5 μ m
- **Flow Rate** : 1.0 mL/min
- **Detection** : PDA 254 nm
- **Temperature** : 35 $^{\circ}$ C
- **Mobile Phase** : Acetonitrile / DW = 60 / 40
- **Samples** : 1) Uracil : 0.1mg/mL
2) Methyl benzoate : 2.2mg/mL
3) Toluene : 8.7mg/mL
4) Naphthalene : 0.9mg/mL

Ghost Peak Removal : Ghost Zero Column


System Connection



Application and Results

| | |
|------------------|---------------------------------|
| Mobile Phase A | Ultra-pure water |
| Mobile Phase B | Acetonitrile |
| Flow Rate | 1.0 mL/min |
| Temperature | 40 °C |
| Detector | 210 nm |
| Injection Volume | 10 µL |
| Sample | Ultra-pure water |
| Column | INNO Column C18 5µm 4.6 * 250mm |

| Time(min) | Mobile Phase A (%) | Mobile Phase B (%) |
|-----------|--------------------|--------------------|
| 0 | 90 | 10 |
| 20 | 10 | 90 |
| 30 | 10 | 90 |
| 30.1 | 90 | 10 |
| 38 | 90 | 10 |




Product Information

| Instrument | Product Name | Column size | Product Number |
|------------|---------------------|-------------|----------------|
| HPLC | Ghost Zero Column | 4.6 x 50mm | GZ04605 |
| UPLC | Ghost Zero Column-U | 2.1 x 30mm | GZ02103-U |
| UPLC | Ghost Zero Column-U | 2.1 x 50mm | GZ02105-U |




HPLC Column : INNO-Core Columns


■ Solid Core Technologies




■ Fast Analysis



■ Retention Time Character of C18 Columns



HPLC Column : INNO-Core Columns



■ Specifications

| Column | Particle Size (µm) | Pore Size (Å) | Surface Area (m ² /g) | C (%) | End Capping | pH Range | USP |
|-------------------------------|--------------------|---------------|----------------------------------|-------|-------------|----------|-----|
| INNO-Core Column C18 | 2.7, 5 | 100 | 115 | 7 | 0 | 2.0~8.5 | L1 |
| INNO-Core Column C18-PE | 2.7, 5 | 100 | 115 | 5 | 0 | 2.0~8.5 | L1 |
| INNO-Core Column C8 | 2.7, 5 | 100 | 115 | 4.5 | 0 | 2.0~8.5 | L7 |
| INNO-Core Column Phenyl-Hexyl | 2.7, 5 | 100 | 115 | 5 | 0 | 2.0~8.5 | L11 |
| INNO-Core Column PFP | 2.7, 5 | 100 | 115 | 5 | 0 | 2.0~8.5 | L43 |

■ Available Particle Size

| Column | 2.7µm | 5µm |
|-------------------------------|-------|-----|
| INNO-Core Column C18 | ● | ● |
| INNO-Core Column C18-PE | ● | ● |
| INNO-Core Column C8 | ● | ● |
| INNO-Core Column Phenyl-Hexyl | ● | ● |
| INNO-Core Column PFP | ● | ● |

■ Available Standard [Horizontal = Inner Diameter (mm) / Vertical = Length (mm)]

| Classify | 2.1 | 3.0 | 4.6 |
|----------|-----|-----|-----|
| 50 | ● | ● | ● |
| 100 | ● | ● | ● |
| 150 | ● | ● | ● |
| 250 | ● | ● | ● |

■ Product Numbering Method

| Particle Size | Product Name (Abbreviation) | Inner Diameter (mm) | Length (mm) |
|-------------------------|---------------------------------------|---|---|
| 2.7µm = 027 5µm = 05 | INNO-Core Column C18 NC | 2.1mm = 021 3.0mm = 030 4.6mm = 046 | 50mm = 05 100mm = 10 150mm = 15 250mm = 25 |
| | INNO-Core Column C18-PE NCPE | | |
| | INNO-Core Column C8 NC8 | | |
| | INNO-Core Column Phenyl-Hexyl NCPH | | |
| | INNO-Core Column PFP NCPFP | | |

| Columns | Particle Size (µm) | Product Name (Abbreviation) | Inner Diameter (mm) | Length (mm) | Product Number |
|--|--------------------|-----------------------------|---------------------|-------------|----------------|
| INNO-Core Column C18 2.7µm 4.6 x 150mm | 027 | NC | 046 | 15 | 027NC04615 |
| INNO-Core Column C18-PE 2.7µm 3.0 x 250mm | 027 | NCPE | 030 | 25 | 027NCPE03025 |
| INNO-Core Column C8 5µm 4.6 x 150mm | 05 | NC8 | 046 | 15 | 05NC804615 |
| INNO-Core Column PFP 2.7µm 2.1 x 100mm | 027 | NCPFP | 021 | 10 | 027NCPFP02110 |



■ Characteristics and parameter

< Aegispak Column Series >

| Column | Particle Size (μm) | Pore Size (Å) | Surface Area (m ² /g) | C(%) | End Capping | pH Range | USP |
|-----------------|--------------------|---------------|----------------------------------|------|-------------|----------|-----|
| Aegispak C18-F | 3,5 | 120 | 320 | 13 | 0 | 1.0~9.0 | L1 |
| Aegispak C18-FA | 3,5 | 120 | 320 | 11 | 0 | 1.0~9.0 | L1 |
| Aegispak C18-L | 3,5 | 100 | 330 | 14 | 0 | 1.0~9.0 | L1 |
| Aegispak C18-LA | 3,5 | 100 | 330 | 13 | 0 | 1.0~9.0 | L1 |
| Aegispak C8 | 3,5 | 100 | 330 | 8 | 0 | 1.0~9.0 | L7 |

< INNO Column Series >

| Column | Particle Size (μm) | Pore Size (Å) | Surface Area (m ² /g) | C(%) | End Capping | pH Range | USP |
|---------------------|--------------------|---------------|----------------------------------|------|-------------|----------|-----|
| INNO Column C18 | 3, 3.5, 5, 10 | 120 | 320 | 18 | 0 | 2.0~8.0 | L1 |
| INNO Column C18-HC | 5 | 100 | 430 | 20 | 0 | 2.0~8.0 | L1 |
| INNO Column C18-SB | 3,5 | 120 | 320 | 17 | 0 | 1.0~11.0 | L1 |
| INNO Column C18-SSB | 5 | 120 | 320 | 18 | 0 | 1.0~11.0 | L1 |
| INNO Column C18-SBA | 3,5 | 120 | 320 | 16 | 0 | 1.0~11.0 | L1 |
| INNO Column C18-NE | 5 | 120 | 320 | 17 | X | 2.0~7.5 | L1 |
| INNO Column C8 | 5 | 120 | 320 | 12 | 0 | 2.0~8.0 | L7 |
| INNO Column C8-HC | 5 | 100 | 430 | 12 | 0 | 2.0~8.0 | L7 |
| INNO Column C8-SB | 5 | 120 | 320 | 11 | 0 | 1.0~11.0 | L7 |
| INNO Column C8-NE | 5 | 120 | 320 | 11 | X | 2.0~7.5 | L7 |
| INNO Column C4 | 5 | 120 | 320 | 8 | 0 | 2.0~8.0 | L26 |
| INNO Column C1 | 5 | 120 | 320 | 3 | X | 2.0~8.0 | L13 |
| INNO Column Silica | 3, 5, 10 | 100, 120 | 320 | - | X | 2.0~8.0 | L3 |
| INNO Column Diol | 5 | 120 | 320 | - | X | 2.0~8.0 | L20 |
| INNO Column NH2 | 3, 5, 10 | 120 | 320 | 6 | X | 2.0~8.0 | L8 |
| INNO Column CN | 5, 10 | 120 | 320 | 8 | 0 | 2.0~8.0 | L10 |
| INNO Column Phenyl | 5 | 120 | 320 | 11 | 0 | 2.0~8.0 | L11 |
| INNO Column SCX | 5 | 120 | 320 | 6 | 0 | 2.0~8.0 | L9 |
| INNO Column SCX-L | 5 | 100 | 320 | 5 | 0 | 2.0~8.0 | L9 |
| INNO Column SAX | 5 | 120 | 320 | 5 | 0 | 2.0~8.0 | L14 |

< INNO-P Column Series >

| Column | Particle Size (μm) | Pore Size (Å) | Surface Area (m ² /g) | C(%) | End Capping | pH Range | USP |
|-------------------|--------------------|---------------|----------------------------------|------|-------------|----------|-----|
| INNO-P Column C18 | 5 | 300 | 110 | 8 | 0 | 1.0~9.0 | L1 |
| INNO-P Column C8 | 5 | 300 | 110 | 5 | 0 | 1.0~9.0 | L7 |
| INNO-P Column C4 | 5 | 300 | 110 | 3 | 0 | 1.0~9.0 | L26 |

< INNO-Mixed Column Series >

| Column | Particle Size (µm) | Pore Size (Å) | Surface Area (m ² /g) | C(%) | End Capping | pH Range | USP |
|---------------------------|--------------------|---------------|----------------------------------|------|-------------|----------|-----|
| INNO Column C18/SCX (1/1) | 5 | 120 | 320 | 12 | 0 | 2.0~8.0 | L## |
| INNO Column C18/SCX (5/1) | 5 | 120 | 320 | 15 | 0 | 2.0~8.0 | L## |
| INNO Column C18/SAX (1/1) | 5 | 120 | 320 | 13 | 0 | 2.0~8.0 | L## |
| INNO Column C18/SAX (5/1) | 5 | 120 | 320 | 16 | 0 | 2.0~8.0 | L## |

< INNO-Core Column Series >

| Column | Particle Size (µm) | Pore Size (Å) | Surface Area (m ² /g) | C(%) | End Capping | pH Range | USP |
|-------------------------------|--------------------|---------------|----------------------------------|------|-------------|----------|-----|
| INNO-Core Column C18 | 2.7, 5 | 100 | 115 | 7 | 0 | 2.0~8.5 | L1 |
| INNO-Core Column C18-PE | 2.7, 5 | 100 | 115 | 5 | 0 | 2.0~8.5 | L1 |
| INNO-Core Column C8 | 2.7, 5 | 100 | 115 | 4.5 | 0 | 2.0~8.5 | L7 |
| INNO-Core Column Phenyl-Hexyl | 2.7, 5 | 100 | 115 | 5 | 0 | 2.0~8.5 | L11 |
| INNO-Core Column PFP | 2.7, 5 | 100 | 115 | 5 | 0 | 2.0~8.5 | L43 |


< Brand Meaning & Model Meaning >

| Brand Name | Brand Meaning | Brand Description |
|-------------------|--------------------------------|--|
| Aegispak | Aegis (Shield) + Pak (Column) | Column capable of analyzing acidic, neutral and basic substances |
| INNO Column | INNO (Innovation) | Separation Revolution Column |
| INNO-P Column | P (Peptide) | Columns of Pore size 300Å |
| INNO-Core Column | C (Core shell) | Columns of Solid Core |
| Ghost Zero Column | Ghost (Peak) + Zero (No exist) | Solvent filter to remove ghosts or unknown peaks |

| Model Name | Model Meaning | Model Description |
|------------|------------------------|---|
| L | Late | Column that elutes later than the F column |
| F | Fast | Column that elute faster than L column |
| FA | Acid | F column treated with strong acid for 6 days |
| LA | Acid | L column treated with strong acid for 6 days |
| SBA | Acid | SB column treated with strong acid for 6 days |
| HC | High Carbon | Columns with maximal binding of C18 or C8 |
| SB | Strong Base | Column with enhanced alkali resistance |
| SSB | Super Strong Base | Column with more enhanced alkali resistance |
| NE | No End capping | Column without end capping |
| PE | Polar Embedded | Column with polar embedded in the C18 |
| PFP | PentaFluoroPhenyl | Column with Pentafluorophenyl group |
| SAX | Strong Anion eXchange | Strong anion exchange column |
| SCX | Strong Cation eXchange | Strong cation exchange column |
| C18/SAX | • | Column with mixed C18 resin and SAX resin |
| C18/SCX | • | Column with mixed C18 resin and SCX resin |




Semi-Prep Column



(HPLC Condition)

- **Column Size** : INNO Column C18
- **Flow Rate** : 1.0 mL/min
4.0 mL/min
9.0 mL/min
- **Detection** : PDA 254 nm
- **Temperature** : 35°C
- **Mobile Phase** : Methanol / DW = 50 / 50
- **Samples** : 1) Uracil : 0.1mg/mL
2) Methyl benzoate : 2.2mg/mL
3) Toluene : 8.7mg/mL
4) Naphthalene : 0.9mg/mL


Allantoin



(HPLC Condition)

- **Column Size** : INNO Column SCX
4.6 x 250m, 5µm
- **Flow Rate** : 1.0 mL/min
- **Detection** : UVD 224nm
- **Temperature** : 40°C
- **Mobile Phase** : Acetonitrile / DW
= 50 / 50
= 70 / 30
= 90 / 10
- **Samples** : Allantoin


Chlorhexidine



(HPLC Condition)

- **Column Size** : INNO Column C18-HC
4.6 x 250m, 5µm
- **Flow Rate** : 1.5 mL/min
- **Detection** : UVD 239nm
- **Temperature** : 40°C
- **Mobile Phase** : Acetonitrile / Buffer=3 / 7
* Buffer : Dissolve 27.6 g of sodium hydrogen phosphate and 10 mL of triethylamine in about 1.5 L of purified water. Add phosphoric acid to pH 3.0 and adjust to 2L.
- **Injection Volume** : 10uL
- **Samples** : 1) Chlorhexidine


Copovidone - USP



(HPLC Condition)

- **Column Size** : Aegispak C18-L 5µm 4.0 x 33mm +
Aegispak C18-L 5µm 4.0 x 250mm
- **Flow Rate** : 1.0 mL/min
- **Detection** : UVD 205nm, 235nm
- **Temperature** : 40°C
- **Mobile Phase** : Acetonitrile /DW = 2 / 23
- **Injection Volume** : 20uL
- **Samples** : 1) 1-Vinyl-2-pyrrolidone
2) Vinyl acetate


Cyclobenzaprine Hydrochloride Tablets - USP



(HPLC Condition)

- **Column Size** : INNO Column C8-SB
4.6 x 250mm, 5 μ m
- **Flow Rate** : 1.0 mL/min
- **Detection** : UVD 226nm
- **Temperature** : 30°C
- **Mobile Phase** : Methanol / Buffer=65 / 35
※ Buffer : 11.4g/L of ammonium acetate in water.
Adjust with ammonium hydroxide to a pH of 7.2.
- **Injection Volume** : 10 μ L
- **Samples** : 1) Cyclobenzaprine Hydrochloride


Hypoalbuminemia



(HPLC Condition)

- **Column Size** : Aegispak C18-L 4.6 X 250mm, 5 μ m
- **Flow Rate** : 0.85 mL/min
- **Detection** : UVD 210 nm
- **Temperature** : 40°C
- **Mobile Phase** : Acetonitrile / pH 2.8 0.2M Sodium Phosphate (mono) = 3/97
- **Samples** : 1) Valine
2) Iso-Leucine
3) Leucine


Lafutidine – JP Method



(HPLC Condition)

- **Column Size** : Aegispak C18-LA 6.0 X 150mm, 5 μ m
- **Flow Rate** : Adjust so that the retention time of lafutidine is about 15 minutes
(About 1.75 mL/min)
- **Detection** : UVD 220 nm
- **Temperature** : 40°C
- **Mobile Phase** : Acetonitrile / Buffer = 15 / 85
- **Buffer** : Dissolve 0.87 g of sodium 1-pentanesulfonate in dilute phosphoric acid (1 in 1000)
- **Samples** : 1) Lafutidine
▶ Theoretical Plate Number=About 10,000

Levocetirizine - EP




(HPLC Condition)

- **Column Size** : INNO Column Silica
4.6 x 250mm, 5 μ m
- **Flow Rate** : 1.0 mL/min
- **Detection** : UVD 230nm
- **Temperature** : 30°C
- **Mobile Phase** : Acetonitrile / DW / 1M H2SO4
= 93 / 6.6 / 0.4
- **Injection Volume** : 20 μ L
- **Samples** : 1) Levocetirizine




Nucleic Acids



(HPLC Condition)

- **Column Size** : Aegispak C18-L 4.6 X 250mm, 5µm
- **Flow Rate** : 1.0 mL/min
- **Detection** : UVD 245 nm
- **Temperature** : 40°C
- **Mobile Phase** : pH 6.4 Phosphate Buffer Solution
- **Samples** : 1) Cytosine
2) Uracil
3) Guanine
4) Thymine
5) Adenine

Organic Acids



(HPLC Condition)

- **Column Size** : Aegispak C18-L 5µm 4.6 X 250mm
- **Flow Rate** : 0.8 mL/min
- **Detection** : PDA 220 nm
- **Temperature** : 30°C
- **Mobile Phase** : Methanol / 10mM Potassium Phosphate Monobasic (pH 2.5, H₃PO₄) = 2 / 98
- **Samples** : 1) Tartaric acid 2) Glycolic acid
3) Malic acid 4) Malonic acid
5) Lactic acid 6) Acetic acid
7) Maleic acid 8) Citric acid
9) Fumaric acid

Pantoprazole Sodium – USP Method




(HPLC Condition)

- **Column Size** : Aegispak C18-L 4.0 X 125mm, 5µm
- **Flow Rate** : 1.0 mL/min
- **Detection** : UVD 290nm and 305nm
- **Temperature** : 40°C
- **Mobile Phase** : Solution (A) = Prepare a Solution of dibasic potassium phosphate(1.74g/L) adjusted with a solution of phosphoric acid(330g/L) to a pH of 7± 0.05
Solution (B) = 100% Acetonitrile
- **Samples** : Pantoprazole
Sodium Related substances test

| Time(minutes) | A | B |
|---------------|-------|-------|
| 0–40 | 80→20 | 20→80 |
| 40–45 | 20→80 | 80→20 |
| 45–55 | 80 | 20 |


Peptide



(HPLC Condition)

- **Column Size** : Aegispak C18-L 4.6 X 250mm, 5µm
- **Flow Rate** : 1.0 mL/min
- **Detection** : PDA 220 nm
- **Temperature** : 40°C
- **Mobile Phase** : Solution A : 0.1 % TFA in Water
Solution B : 0.1 % TFA in (Acetonitrile / Methanol = 50/50)
- **Samples** : 1) Pentapeptide 2) Caffeine
3) 15n-labeled amino acid
4) 10n-labeled amino acid
5) Coumaric acid
6) 15n-labeled amino acid + Coumaric acid


Pranlukast Hydrate - JP Method



(HPLC Condition)

- **Column Size** : Aegispak C8 6.0 X 150mm, 5µm
- **Flow Rate** : Adjust so that the retention time of pranlukast is about 10 minutes
- **Detection** : UVD 260 nm
- **Temperature** : 25°C
- **Mobile Phase** : A mixture of 0.02 mol/L potassium dihydrogen phosphate TS, acetonitrile and methanol(5:5:1)
- **Samples** : Pranlukast hydrate
Assay test


Preservatives



(HPLC Condition)

- **Column Size** : Aegispak C18-F 5µm 4.6 X 150mm
- **Mobile Phase** : A/B=35/65 A : Acetonitrile
B : 50mM Ammonium Acetate Buffer Soln.(pH 4.6)
- **Flow Rate** : 1.0ml/min
- **Wavelength** : PDA 235 nm
- **Oven Temp** : 40°C
- **Samples** : 1. Benzoic acid 2. Sorbic acid
3. Dehydroacetic acid 4. Methyl Parabene
5. Ethyl Parabene 6. Isopropyl Parabene
7. Propyl Parabene 8. Isobutyl Parabene
9. Butyl Parabene

Protamine Sulfate - USP Method




(HPLC Condition)

- **Column Size** : INNO-P Column C18 4.6 X 250mm, 5µm
- **Flow Rate** : 1.0 mL/min
- **Detection** : UVD 214 nm
- **Temperature** : 40°C
- **Mobile Phase** : Solution A : 0.1 M Sodium dihydrogen phosphate solution (adjust with pH 1.8 phosphoric acid)
Solution B : Solution A 467.5 mL + Acetonitrile 32.5 mL

| Time(min) | Solution A | Solution B |
|-----------|------------|------------|
| 0 | 85 | 15 |
| 15 | 55 | 45 |
| 25 | 55 | 45 |
| 30 | 85 | 15 |

Rebamipide - JP Method



(HPLC Condition)

- **Column Size** : INNO Column C18-SB 4.6 X 150mm, 5µm
- **Flow Rate** : Adjust so that the retention time of rebamipide is about 20 minutes (about 1.1 mL/min)
- **Detection** : UVD 222 nm
- **Temperature** : 25°C
- **Mobile Phase** : To 300 mL of phosphate buffer solution (pH 6.2) add 750 mL of water. To 830 mL of this solution add 170 mL of acetonitrile
- **Samples** : Rebamipide m-chloro isomer test
▶ Theoretical Plate Number = About 14,000



Rosuvastatin calcium - EP Method




(HPLC Condition)

- **Column Size** : INNO Column C18 3.0 X 150mm, 3µm
- **Flow Rate** : 0.75 mL/min
- **Detection** : UVD 242 nm
- **Temperature** : 40°C
- **Mobile Phase** : Solution A: 1 % Trifluoroacetic acid : Acetonitrile : water (1:29:70)
Solution B: 1 % Trifluoroacetic acid : water : Acetonitrile (1:24:75)
- **Samples** : Rosuvastatin calcium / Related substances test

| Time(min) | Solution A | Solution B |
|-----------|------------|------------|
| 0-30 | 100 | 0 |
| 30-50 | 100→60 | 0→40 |
| 50-60 | 60→0 | 40→100 |
| 60-70 | 0 | 100 |


Tocopherol calcium succinate - KP Method



(HPLC Condition)

- **Column Size** : Aegispak C18-L 4.6 X 150mm, 5µm
- **Flow Rate** : Adjust so that retention time of tocopherol succinate is about 8 minutes (about 0.75 mL/min)
- **Detection** : UVD 284 nm
- **Temperature** : 25°C
- **Mobile Phase** : Methanol : water : acetic acid(100) (97:2:1)
- **Samples** : 1) Tocopherol succinate
2) Tocopherol


Vitamin C (Ascorbic acid)



(HPLC Condition)

- **Column Size** : INNO Column C18-HC 4.6 X 250mm, 5µm
INNO Column C18-NE 4.6 X 250mm, 5µm
INNO Column NH2 4.6 X 250mm, 5µm
INNO Column C18-HC 4.0 X 250mm, 5µm
- **Flow Rate** : 1.0, 1.0, 1.0, 0.75 mL/min
- **Detection** : UVD 254 nm
- **Temperature** : 25°C
- **Mobile Phase** : Acetonitrile / 0.05M KH₂PO₄ = 40 / 60
0.05 M KH₂PO₄ : Dissolve 6.8 g of potassium dihydrogen phosphate in water to make 1000 mL.
- **Samples** : 1) Ascorbic acid

Vitamins – Water soluble



(HPLC Condition)

- **Column Size** : INNO Column C18 5µm 4.6 X 250mm
- **Flow Rate** : 0.8 mL/min
- **Detection** : PDA 210 nm
- **Temperature** : 35°C
- **Mobile Phase** : A) 10mM KH₂PO₄ Buffer(pH 3.0, H₃PO₄) Acetonitrile = 98/2
B) 10mM KH₂PO₄ (pH 3.0, H₃PO₄) Acetonitrile = 40/60
- **Samples** : 1) Thiamine 2) Ascorbic acid 3) Nicotinic acid
4) Pyridoxine 5) Pyridoxal 6) Nicotinamide
7) Pantothenic acid 8) Hydroxocobalamin
9) Cyanocobalamin 10) Folic acid
11) Riboflavin 12) Biotin

Product List

Available particle size

| Available Particle Size | | 2.7µm | 3µm | 3.5µm | 5µm | 10µm(100 Å) | 10µm(120 Å) |
|-------------------------------|-------|-------|-----|-------|-----|-------------|-------------|
| Aegispak C18-F | AF | | • | | • | | |
| Aegispak C18-FA | AFA | | • | | • | | |
| Aegispak C18-L | AL | | • | | • | | |
| Aegispak C18-LA | ALA | | • | | • | | |
| Aegispak C8 | A8 | | • | | • | | |
| INNO Column C18 | N | | • | • | • | • | • |
| INNO Column C18-HC | NHC | | | | • | | |
| INNO Column C18-SB | NSB | | • | | • | | |
| INNO Column C18-SSB | NSSB | | | | • | | |
| INNO Column C18-SBA | NSBA | | | | • | | |
| INNO Column C18-NE | NE | | • | | • | | |
| INNO Column C8 | N8 | | | | • | | |
| INNO Column C8-HC | N8HC | | | | • | | |
| INNO Column C8-SB | N8SB | | | | • | | |
| INNO Column C8-NE | N8E | | | | • | | |
| INNO Column C4 | N4 | | | | • | | |
| INNO Column C1 | N1 | | | | • | | |
| INNO Column Phenyl | PH | | | | • | | |
| INNO Column Silica | S | | • | | • | • | • |
| INNO Column Diol | DL | | | | • | | |
| INNO Column NH2 | NH | | • | | • | | • |
| INNO Column CN | CN | | | | • | | • |
| INNO Column SCX | SCX | | | | • | | |
| INNO Column SCX-L | SCX-L | | | | • | | |
| INNO Column SAX | SAX | | | | • | | |
| INNO Column C18/SCX (1/1) | N1SCX | | | | • | | |
| INNO Column C18/SCX (5/1) | N5SCX | | | | • | | |
| INNO Column C18/SAX (1/1) | N1SAX | | | | • | | |
| INNO Column C18/SAX (5/1) | N5SAX | | | | • | | |
| INNO-P Column C18 | NP | | | | • | | |
| INNO-P Column C8 | NP8 | | | | • | | |
| INNO-P Column C4 | NP4 | | | | • | | |
| INNO-Core Column C18 | NC | • | | | • | | |
| INNO-Core Column C18-PE | NCPE | • | | | • | | |
| INNO-Core Column C8 | NC8 | • | | | • | | |
| INNO-Core Column Phenyl-Hexyl | NCPH | • | | | • | | |
| INNO-Core Column PFP | NCPFP | • | | | • | | |





■ Available Standard (Horizontal = Inner Diameter (mm)/Vertical = Length (mm))

| Classify | 1.5 | 2.0 | 2.1 | 3.0 | 3.9 | 4.0 | 4.6 | 6.0 | 8.0 | 10 | 20 |
|----------|-----|-----|-----|-----|-----|-----|-----|-----|-----|----|----|
| 30 | • | • | | • | • | • | • | | | | |
| 33 | • | • | | • | • | • | • | | | | |
| 35 | • | • | | • | • | • | • | | | | |
| 40 | • | • | | • | • | • | • | | | | |
| 45 | • | • | | • | • | • | • | | | | |
| 50 | • | • | • | • | • | • | • | | | | |
| 55 | • | • | | • | • | • | • | | | | |
| 60 | • | • | | • | • | • | • | | | | |
| 75 | • | • | | • | • | • | • | | | | |
| 80 | • | • | | • | • | • | • | | | | |
| 100 | • | • | • | • | • | • | • | • | • | | |
| 110 | • | • | | • | • | • | • | | | | |
| 120 | • | • | | • | • | • | • | | | | |
| 125 | • | • | | • | • | • | • | | | | |
| 150 | • | • | • | • | • | • | • | • | • | | |
| 200 | • | • | | • | • | • | • | • | • | | |
| 220 | • | • | | • | • | • | • | • | • | | |
| 250 | • | • | • | • | • | • | • | • | • | • | • |
| 300 | | | | • | • | • | • | • | • | | |

※ In the available specification, the length can be supplied as you wish.
 ※ Inner diameter 2,1mm is only for INNO-Core Column.

Product numbering method

| Particle size (um) | Product Name (Abbreviation) | Inner Diameter (mm) | Length (cm) |
|---|---|--|---|
| 03 | AL | 015 | 03 |
| 2.7um = 027 3um = 03 3.5um = 035 5um = 05 10um = 10 (P/S = 120Å) 10um=010 (P/S = 100Å) | Aegispak C18-F = AF Aegispak C18-FA= AFA Aegispak C18-L = AL Aegispak C18-LA = ALA Aegispak C8 = A8 INNO Column C18 = N INNO Column C18-HC = NHC INNO Column C18-SB = NSB INNO Column C18-SSB = NSSB INNO Column C18-SBA = NSBA INNO Column C18-NE = NE INNO Column C8 = N8 INNO Column C8-HC = N8HC INNO Column C8-SB = N8SB INNO Column C8-NE = N8E INNO Column C4 = N4 INNO Column C1 = N1 INNO Column Silica = S INNO Column Diol = DL INNO Column CN = CN INNO Column NH2 = NH INNO Column Phenyl = PH INNO Column SCX = SCX INNO Column SCX-L = SCXL INNO Column SAX = SAX INNO Column C18/SCX (1/1) = N1SCX INNO Column C18/SCX (5/1) = N5SCX INNO Column C18/SAX (1/1) = N1SAX INNO Column C18/SAX (5/1) = N5SAX INNO-P Column C18 = NP INNO-P Column C8 = NP8 INNO-P Column C4 = NP4 INNO-Core Column C18 = NC INNO-Core Column C18-PE = NCPE INNO-Core Column C8 = NC8 INNO-Core Column Phenyl-Hexyl = NCPH INNO-Core Column PFP = NCPFP | 1.5mm = 015 2.0mm = 020 3.0mm = 030 3.9mm = 039 4.0mm = 040 4.6mm = 046 6.0mm = 060 8.0mm = 080 10.0mm = 100 20.0mm = 200 30.0mm = 300 | 30mm = 03 33mm = 033 35mm = 035 40mm = 04 45mm = 045 50mm = 05 55mm = 055 60mm = 06 70mm = 07 75mm = 075 80mm = 08 100mm = 10 110mm = 11 120mm = 12 125mm = 125 150mm = 15 200mm = 20 220mm = 22 250mm = 25 300mm = 30 |

Examples

| Columns(um) | Particle Size (um) | Product Name (Abbreviation) | Inner Diameter(mm) | Length (cm) | Product NO |
|--|--------------------|-----------------------------|--------------------|-------------|-------------|
| Aegispak C18-L 5um 3.9 x 150mm | 05 | AL | 039 | 15 | 05AL03915 |
| Aegispak C8 5um 4.6 x 250mm | 05 | A8 | 046 | 25 | 05A804625 |
| INNO Column C18-SB 5um 4.6 x 125mm | 05 | NSB | 046 | 12.5 | 05NSB046125 |
| INNO Column CN 5um 4.0 x 250mm | 05 | CN | 040 | 25 | 05CN04025 |
| INNO Column NH2 5um 4.6 x 250mm | 05 | NH | 046 | 25 | 05NH04625 |
| INNO Column SCX 5um 4.6 x 250mm | 05 | SCX | 046 | 25 | 05SCX04625 |
| INNO-P Column C4 5um 4.6 x 250mm | 05 | NP4 | 046 | 25 | 05NP404625 |
| INNO-Core Column C18 2.7um 4.6 x 150mm | 027 | NC | 046 | 15 | 027NC04615 |



■ USP L Code

| L Code | | Columns |
|--------|--|-------------------------------|
| L1 | Octadecyl silane chemically bonded to porous or non-porous silica or ceramic micro-particles, 1.5 to 10 µm in diameter, or a monolithic rod | Aegispak C18-F |
| | | Aegispak C18-FA |
| | | Aegispak C18-L |
| | | Aegispak C18-LA |
| | | INNO Column C18 |
| | | INNO Column C18-HC |
| | | INNO Column C18-SB |
| | | INNO Column C18-SSB |
| | | INNO Column C18-SBA |
| | | INNO Column C18-NE |
| | | INNO-P Column C18 |
| | | INNO-Core Column C18 |
| | | INNO-Core Column C18-PE |
| L3 | Porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod. | Aegispak Silica |
| | | INNO Column Silica |
| L7 | Octylsilane chemically bonded to totally or superficially porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod. | Aegispak C8 |
| | | INNO Column C8 |
| | | INNO Column C8-HC |
| | | INNO Column C8-SB |
| | | INNO Column C8-NE |
| | | INNO-P Column C8 |
| | | INNO-Core Column C8 |
| L8 | An essentially monomolecular layer of aminopropylsilane chemically bonded to totally porous silica gel support, 1.5 to 10 µm in diameter, or a monolithic silica rod. | INNO Column NH2 |
| L9 | Irregular or spherical, totally porous silica gel having a chemically bonded, strongly acidic cation-exchange coating, 3 to 10 µm in diameter | INNO Column SCX |
| | | INNO Column SCX-L |
| L10 | Nitrile groups chemically bonded to porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod.. | INNO Column CN |
| L11 | Phenyl groups chemically bonded to porous silica particles 1.5 to 10 µm in diameter, or a monolithic silica rod. | INNO Column Phenyl |
| | | INNO-Core Column Phenyl-Hexyl |
| L13 | Trimethylsilane groups chemically bonded to porous silica or hybrid particles, 1.5 to 10 µm in diameter | INNO Column C1 |
| L14 | Silica gel having a chemically bonded, strongly basic quaternary ammonium anionexchange coating, 5 to 10 µm in diameter. | INNO Column SAX |
| L20 | Dihydroxypropane groups chemically bonded to porous silica or hybrid particles, 1.5 to 10 µm in diameter. | INNO Column Diol |
| L26 | Butyl silane chemically bonded to totally porous silica particles, 1.5 to 10 µm in diameter | INNO Column C4 |
| | | INNO-P Column C4 |
| L43 | Pentafluorophenyl groups chemically bonded to porous or superficially porous silica particles by a propyl spacer, 1.5 to 10 µm in diameter. | INNO-Core Column PFP |
| L## | Octadecyl silane chemically bonded to porous silica particles, 1.5 to 10 µm in diameter and porous silica gel having a chemically bonded, strongly acidic cation-exchange coating, 3 to 10 µm in diameter. The two resins are mixed 1/1 or 5/1 | INNO Column C18/SCX |
| L## | Octadecyl silane chemically bonded to porous silica particles, 1.5 to 10 µm in diameter and porous silica gel having a chemically bonded, strongly basic quaternary ammonium anion exchange coating, 3 to 10 µm in diameter. The two resins are mixed 1/1 or 5/1 | INNO Column C18/SAX |
| L## | This column is a unit that removes unknown or ghost peaks from the mobile phase. | Ghost Zero Column |



■ Correspondence Table

| Youngjin Biochrom | INNO-P | Aegispak C18-FA | Aegispak C18-F ----- Aegispak C18-LA | Aegispak C18-L ----- INNO Column C18-SSB | INNO Column C18 | INNO Column C18-SB / C18-SSB / C18-SBA | INNO Column C18-HC | INNO Column C18-NE | INNO-Core Column C18 |
|-------------------|------------------------------------|----------------------------|--|---|--|--|--|--|----------------------------|
| Agilent | ZORBAX 300SB-C18 | ZORBAX SB-C18 | Zorbax Eclipse C18 | Zorbax Eclipse Plus C18 | | | | | Poroshell 120 C18 |
| AMT | | | | | | | | | Halo C18 |
| EKA | | | | Kromasil Eternity C18 | Kromasil C18 | | | | Kromasil EternityShell C18 |
| GL Science | Inertsil ODS-SP | InertSustain Swift C18 | Inertsil ODS-2 | Inertsil ODS-3 Inertsil ODS-4 InertSustain C18 | | | Inertsil ODS-P | Inertsil ODS-EP | InertCore Plus C18 |
| Merck | LiChrospher® WP 300 RP-18 | | | | Purospher (STAR) RP-18 endcapped | | LiChrospher 100 RP-18 endcapped LiChrospher PAH Purospher RP-18 HC | LiChrosorb RP-18 | Ascentis® Express C18 |
| Nacalai | Cosmosil C18 AR-300 | Cosmosil C18-PAQ | Cosmosil C18-EB | Cosmosil C18 MS-II | Cosmosil C18 AR-II | | | | COSMOCORE C18 |
| Nagel | Nucleosil 300 C18 | | | Nucleosil C18 | Nucleodur C18 Gravity Nucleodur C18 Htec | | Nucleodur C18 ISIS Nucleosil C18 HD | Nucleosil C18 Nautilus Nucleodur C18 Pyramid | Nucleoshell™ RP-C18 |
| Nomura | | Develosil ODS-P | Develosil ODS-UG Develosil ODS-HG | Develosil ODS-MG Develosil XG C18LC | Develosil ODS-SR Develosil XG C18M | | Develosil ODS | | |
| OSAKA Soda | | | Capcellpak C18 UG120 | Capcellpak C18 MG Capcellpak C18 AQ | | | | | Capcellcore C18 |
| Phenomenex | | | | Gemini C18 Gemini NX C18 | Luna C18 (2) Onyx C18 (HD) | EVO C18 | | | Kinetex C18 |
| Supelco | Discovery BIO Wide Pore C18 | | SUPELCOSIL™ LC-18 | SUPELCOSIL™ LC-18-DB | Discovery C18 | | Discovery HS C18 | | |
| Thermo | | | Hypersil ODS | Hypersil BDS Hypersil Gold C18 | Acclaim 120 C18 | | | | Accucore RP-MS |
| Vydac | Vydac 201TP C18 Vydac 218MS C18 | | | | | | | | |
| Waters | | X-Terra XSelect HSS C18 SB | XTerra Shield RP18 Spherisorb ODS-2 | XTerra MS C18 Sunfire C18 X-Bridge C18 XSelect CSH C18 | Athlantis C18 Symmetry (Shield) C18 XBridge BEH C18 uBondapak C18 | X-Bridge C18 | | Spherisorb ODS-1 | Cortecs Solid Core Column |
| YMC | | | YMCpak Hydrosphere C18 | YMCpak Pro C18 | YMCpak ODS-AQ | YMCpak Triart C18 | YMCpak Triart C18 | | Meteoric Core Column C18 |

The corresponding columns may behave differently depending on the analysis conditions and the analyte.





YJ BIOCHROM
<http://www.yjbiochrom.com>

1301, Dunchon-daero 388, Jungwon-gu,
Seongnam-si, Gyeonggi-do, Korea
TEL : +82-31-736-4103
FAX : +82-31-736-4104
<http://www.yjbiochrom.com>

